

### Safe Harbor Statement

This presentation contains "forward-looking statements" pursuant to the safe harbor provisions of the Private Securities Litigation Reform Act of 1995. Forward-looking statements include, without limitation, any statement that may predict, forecast, indicate, or imply future results, performance or achievements, and may contain the words "estimate," "intend," "target," "will," "is likely," "would," "may," or, in each case, their negative, or words or expressions of similar meaning. These forward-looking statements are found at various places throughout this presentation and include information concerning possible or assumed future results of our operations; business strategies; future cash flows; financing plans; plans and objectives of management; any other statements regarding future operations, future cash needs, business plans and future financial results, and any other statements that are not historical facts. Unless otherwise indicated, the terms "CytoSorbents," "Company," "we," "us" and "our" refer to CytoSorbents Corporation. Any or all of the forwardlooking statements included in this presentation are not guarantees of future performance and may turn out to be inaccurate. These forward-looking statements represent our intentions, plans, expectations, assumptions and beliefs about future events and are subject to risks, uncertainties and other factors. Many of those factors are outside of our control and could cause actual results to differ materially from the results expressed or implied by those forward-looking statements. Although these expectations may change, we are under no obligation to inform you if they do. Actual events or results may differ materially from those contained in the forward-looking statements. The following factors, among others, could cause our actual results to differ materially from those described in a forward-looking statement: our history of losses; potential fluctuations in our quarterly and annual results; competition, inability to achieve regulatory approval for our device, technology systems beyond our control and technology-related defects that could affect the companies' products or reputation; risks related to adverse business conditions; our dependence on key employees; competition for qualified personnel; the possible unavailability of financing as and if needed; and risks related to protecting our intellectual property rights or potential infringement of the intellectual property rights of third parties. This list is intended to identify only certain of the principal factors that could cause actual results to differ from those discussed in the forward-looking statements. In light of these risks, uncertainties and assumptions, the events described in the forwardlooking statements might not occur or might occur to a different extent or at a different time than we have described. You are cautioned not to place undue reliance on these forward-looking statements, which speak only as of the date of the applicable presentation. You are referred to a discussion of important risk factors detailed in the Company's Form 10-K filed with the Securities and Exchange Commission on March 5, 2020 and other reports and documents filed from time to time by us, which are available online at www.sec.gov.

## **Opening Remarks**

Dr. Phillip Chan, MD, PhD Chief Executive Officer CytoSorbents Corporation

## **CytoSorbents** At a Glance (NASDAQ: CTSO)

- CytoSorbents is a U.S. NASDAQ-traded medical device company that specializes in treating life-threatening conditions with its blood purification technology
- CytoSorb® is E.U. approved, manufactured in the U.S. by CytoSorbents, and commercialized in 65 countries as an extracorporeal cytokine adsorber to help treat hyperinflammatory conditions where cytokines are elevated (e.g. "cytokine storm") with more than 88,000 cumulative treatments to date
- CytoSorb is also E.U. approved to remove ticagrelor (Brilinta®) or rivaroxaban (Xarelto®) in cardiac surgery, bilirubin (liver disease) and myoglobin (trauma)
- CytoSorb is not yet FDA-approved but on a dual path for U.S. approval
  - FDA Breakthrough Designation to remove ticagrelor during CPB in urgent & emergent cardiothoracic surgery
  - U.S. REFRESH 2-AKI Trial Pivotal study at 25 U.S. centers using CytoSorb intraoperatively to reduce risk of post-op AKI
- Received U.S. FDA Emergency Use Authorization for use in critically-ill adult COVID-19+ patients with respiratory failure and has been used ~1,000 COVID-19 patients in 20+ countries, including the U.S.
- 156 employees with international footprint across two wholly-owned subsidiaries
  - CytoSorbents Medical, Inc: Headquarters New Jersey, USA (ISO 13485 certified manufacturing, R&D, Management)
  - CytoSorbents Europe GmbH: International sales office Berlin, Germany (Sales and Marketing)
- Strong government support with ~\$33M in grants, contracts, other non-dilutive funds




## CytoSorb is "Plug and Play" Compatible

#### Compatible with Existing Blood Pump Infrastructure In Hospitals Today


#### **Dialysis or CRRT**

(Continuous Renal Replacement Therapy)




#### Hemoperfusion

(Standalone Treatment)




#### **ECMO**

(Extracorporeal Membrane Oxygenation)





#### **CPB**

(Cardiopulmonary Bypass)



### **CytoSorb** Commercialization Focus





#### **World Class Partners**



(Critical Care)



(Cardiac Surgery)




(India – Cardiac / Critical Care)



## CytoSorb Adoption Continues to Grow

#### Product Sales(ttm) and Blended Product Gross Margin Growth




## Growth Driven By Many Macro Trends in Healthcare




































## Introduction of Agenda & Speakers

Efthymios N. Deliargyris, MD, FACC, FESC, FSCAI Chief Medical Officer CytoSorbents Corporation



## CytoSorb Drug Adsorption: Mechanism of Action



- favors hydrophobic binding
- pore size = selective access to ≤ 60 kDa
- dependent on drug concentration
- dependent on time of blood exposure


### **Drug:**

- Hydrophobic molecular functional groups (aromatic, alkyl, etc.) increase binding affinity
- Free fraction vs. protein-bound: free fraction adsorbed and shifting free-bound ratio
- Active metabolites: generation rate and compartmentalization factor into removal rates
- $T_{1/2}$ : very short  $T_{1/2}$  limits availability for adsorption vs. longer  $T_{1/2}$
- Volume of distribution (V<sub>D</sub>): more adsorption with smaller V<sub>D</sub>
- Chronicity of dosing: higher adsorption early before steady-state V<sub>D</sub>




### CytoSorb Integrates Easily Into Cardiopulmonary Bypass

- CytoSorb installs within minutes and is placed in a parallel circuit: post-pump, back to the venous reservoir
- High blood flow, low resistance up to 700 mL/min
- Fully-compatible with heparin anti-coagulation
- Used safely in thousands of cardiopulmonary bypass procedures to date





### **Antithrombotic Removal - Overview**





## Today's Agenda & Speakers

Speaker	Title	Time
Phillip Chan	Welcome & Opening Remarks	5 min
Makis Deliargyris	Introduction of Agenda & Speakers	5 min
Robert Storey	Ticagrelor and CytoSorb	15 min
Michael Schmoeckel	Intraoperative removal of Ticagrelor and Rivaroxaban during Emergency Cardiac Operations	15 min
Michael Gibson	NOACs and CytoSorb	15 min
Makis Deliargyris	Closing Remarks – Size of the opportunity	5 min
All	Q & A Session	30 min







## **Ticagrelor and CytoSorb**


Professor Robert Storey, BSc, BM, DM, FRCP, FESC

Professor of Cardiology and Cardiovascular Disease Theme Lead, Department of Infection, Immunity and Cardiovascular Disease, University of Sheffield

and


Academic Director and Honorary Consultant Cardiologist, Cardiology and Cardiothoracic Surgery Directorate, Sheffield Teaching Hospitals NHS Foundation Trust, Sheffield, United Kingdom

# Ticagrelor: the first oral reversibly-binding P2Y<sub>12</sub> receptor antagonist belonging to the class CPTP (cyclo-pentyl-triazolo-pyrimidine)




van Giezen JJJ, Humphries RG. Semin Thromb Hemost. 2005;31:195-204.

#### **ONSET/OFFSET Study: inhibition of platelet aggregation**




## **PLATO PLATELET: VerifyNow P2Y<sub>12</sub> Assay**

**Comparing Maintenance Therapy with Clopidogrel vs Ticagrelor in ACS** 



\*\*\*P<0.0001; PRU = Platelet reaction units. Storey RF, et al. J Am Coll Cardiol 2010

### **PLATO Secondary efficacy endpoints**



Wallentin L et al. N Eng J Med. 2009


## UK networks for ACS and revascularisation

## **South Yorkshire Cardiothoracic Centre**

Provides a PCI and CABG surgery service including 24/7 primary PCI to the South Yorkshire and North Derbyshire regions of England – a population of 1.8 million people




## Sheffield observational study 10,793 consecutive invasively-managed ACS patients



Iqbal J et al. Presented at AHA 2014; Gosling R et al. Platelets 2017 online

## Sheffield observational study Adjusted definite stent thrombosis rates






# BCIS National Audit Adult Interventional Procedures

1st April 2018 to 31st March 2019






## Clinical Syndrome








## Ticagrelor Use by Indication for PCI





#### Isolated CABG surgery: elective vs. urgent




Mean time waiting for urgent CABG in 2017/18 was 10 days in the UK (no change from 2016/17)


Proportion treated within 7 days = 34%

Only 4 hospitals managed to treat >50% within 7 days


# PLATO: chest tube drainage according to time after drug intake



# PLATO: fall in hemoglobin >50g/L according to time after drug intake



# Offset of ticagrelor's effects in ACS patients



\*\*\*



#### Expert position paper on the management of antiplatelet therapy in patients undergoing coronary artery bypass graft surgery

Miguel Sousa-Uva<sup>1,2</sup>, Robert Storey<sup>3</sup>, Kurt Huber<sup>4</sup>, Volkmar Falk<sup>5</sup>, Adeline Leite-Moreira<sup>6,7</sup>, Julien Amour<sup>8</sup>, Nawwar Al- Attar<sup>9</sup>, Raimondo Ascione<sup>10</sup>, David Taggart<sup>11</sup>, and Jean-Philippe Collet<sup>8\*</sup>, on behalf of ESC Working Group on Cardiovascular Surgery and ESC Working Group on Thrombosis

Table 3 Proposed strategies for discontinuation of P2Y<sub>12</sub> inhibitors prior to coronary artery bypass grafting surgery


	Bleeding risk							
Thrombotic risk	High <sup>a</sup>		Low					
	High <sup>b</sup>	Early Heart Team Consultation	Early Heart Team Consultation					
	ACS or recent stent PCI	Ticagrelor/clopidogrel: stop 5 days before and bridge for 4 days. Prasugrel: stop 7 days and bridge for 5 days	Ticagrelor/clopidogrel: stop 3 days before and bridge for 2 days. Prasugrel: stop 5 days before and bridge for 3 days					
	Low	Early Heart Team Consultation Clopidogrel/ ticagrelor: stop 5 days before. Prasugrel: stop 7 days prior to CABG	Clopidogrel/ticagrelor: stop 5 days before or less if indicated by platelet function test. Prasugrel: stop 7 days before or less if indicated by platelet function test.					

<sup>&</sup>lt;sup>a</sup>Examples of high-bleeding risk: renal or hepatic insufficiency, advanced age, anaemia, small body surface area, cardiac failure, and redoes operation.

<sup>&</sup>lt;sup>b</sup>Examples of high-thrombotic risk: haemodynamic instability, ongoing ischaemia, complex coronary anatomy, stenting < 1 month for BMS, and < 6 months for DES. CABG, coronary atery bypass grafting.

## **Ticagrelor Removal From Human Blood**

George O. Angheloiu, MD,<sup>a,b,c</sup> Gabriel B. Gugiu, PнD,<sup>d</sup> Cristian Ruse, PнD,<sup>e</sup> Rishikesh Pandey, PнD,<sup>a</sup> Ramachandra R. Dasari, PнD,<sup>a</sup> Carl Whatling, PнD<sup>f</sup>



Angheloiu GO et al. JACC Basic Trans Sci 2017



# Ticagrelor removal by CytoSorb® is associated with reduced morbidity in patients who require emergent or urgent cardiac surgery:

An economic model with implications for hospital resource utilisation in the UK

M. Javanbakht<sup>1</sup>, K. Rahimi<sup>2</sup>, F. Degener<sup>3</sup>, D. Adam<sup>3</sup>, F. Preissing<sup>3</sup>, J. Scheier<sup>4</sup>, SF. Cook<sup>5</sup>, **E. Mortensen**<sup>6\*</sup>

1. Optimax Access UK Ltd, Market Access Consultancy, Southampton, United Kingdom; 2. The George Institute for Global Health, University of Oxford, Oxford, United Kingdom; 3. Reimbursement & Health Economics, CytoSorbents Europe GmbH, Berlin, Germany; 4. Medical Affairs, CytoSorbents Europe GmbH, Berlin, Germany; 5. CERobs Consulting LLC, Chapel Hill, United States of America; 6. Medical Affairs, CytoSorbents Corporation, Monmouth Junction, United States of America \*Presenting author

#### Methods

- A de novo decision analytic model was developed to estimate resource utilisation in each strategy (CytoSorb vs. usual care) over a 30-day time horizon
- Primary clinical inputs were those that might have significant impact on hospital care resource utilisation, including:
  - bleeding complications
  - re-thoracotomies
  - number of transfused units of red blood cells (RBC) and platelets
  - hospital/intensive care unit length of stay and total operating time
  - incidence of myocardial infarction while waiting for physiologic clearance of ticagrelor before an urgent cardiac surgery
- A wide range of parametric and structural sensitivity analyses were performed to explore the uncertainty surrounding the model results

#### Results: Emergent Cardiac Surgery (Cohort 1)

- CytoSorb use resulted in fewer blood product transfusions, fewer re-thoracotomies, shorter operation time, and shorter ICU/hospital length of stay
- Patients treated with CytoSorb incurred lower total cost (-£3,982) and had better health-related quality of life (+0.00125 quality-adjusted life years (QALYs))
- CytoSorb remained dominant in all sensitivity analyses

Outcome (30-day time horizon;	Without	With	Δ
Comparator: No physiologic clearance)	CytoSorb	CytoSorb	incremental
RBCs transfusions (units) per patient	0.91	0.44	-0.47
Platelet transfusions (units) per patient	1.55	0.38	-1.16
Operation time (in minutes)	353	288	-65
Re-thoracotomy rate	36%	0%	-100%
ICU length of stay (days)	3	2	-1
Hospital length of stay (days)	14	11	-3

Javanbakht M et al. Presented at EACTS 2019 meeting, Lisbon, Portugal, October 2019

#### Results: Urgent Cardiac Surgery (Cohort 2)


- In urgent CABG, CytoSorb was less costly (-£55) and more effective when compared to waiting for
   5 days to allow for physiological washout of ticagrelor to reduce bleeding risk
- In all 3 comparators, waiting alone, waiting plus short acting antiplatelet agent, and waiting plus low molecular weight heparin:
  - As expected, transfusion of blood and platelet are similar with or without CytoSorb treatment, as bleeding risk is reduced in both cases
  - Hospital length of stay was reduced by 5 days as surgery could proceed earlier

Outcome per patient (Deterministic)	Without CytoSorb	With CytoSorb	Δ Incre- mental
Transfusions of red blood cells and platelets and (units) per patient	0.82	0.82	0
Hospital length of stay (days)	16	11	-5


#### <u>Ti</u>cagrelor Cyto<u>Sorb</u> Hemoadsorption (TISORB) study

UK prospective, multi-center study in patients undergoing cardiac surgery <48 hours since last dose of ticagrelor






#### TISORB patient journey: Surgery ≤ 48 hrs after last ticagrelor dose



#### **TISORB** primary effectiveness endpoint



#### **Summary**

- Acute coronary syndromes (ACS) have become the dominant reason for revascularization over the last decade
- Dual antiplatelet therapy with aspirin and ticagrelor is first-line therapy for ACS patients and has dramatically cut stent thrombosis rates in PCI patients
- Level of platelet P2Y<sub>12</sub> inhibition has a critical effect on surgical blood loss and the risks of urgent CABG surgery
- Ticagrelor has the advantage over irreversible P2Y<sub>12</sub> inhibitors of reversibility and can be removed from the blood by the CytoSorb system
- CytoSorb has the potential to transform the safety, timeliness, simplicity and costeffectiveness of CABG surgery in the ACS population



# Intraoperative removal of Ticagrelor and Rivaroxaban during Emergency Cardiac Operations




Prof. Dr. med. Michael Schmoeckel
Head
Dept. of Cardiac Surgery
AK St. Georg, Hamburg, Germany



#### **Perioperative management of NOACs**



- Apixaban (Eliquis) discontinue apixaban 24 to 48 hours prior to surgery depending on the bleeding risk.
- Dabigatran (Pradaxa) high bleeding risk procedures or surgeries:
   to be discontinued 48–72 hours before.
   in renal impairment: 72–96 hours before.
- Rivaroxaban (Xarelto) to be discontinued 48 hours prior to high bleeding risk procedures.
- Edoxaban (Savaysa) in high bleeding risk procedures, edoxaban should be discontinued 72 hours before.



between July 2014 and June 2016. All patients presented for surgery while on NOAC therapy: 37 received rivaroxaban (45.7%), 35 apixaban (43.2%), and 9 dabigatran (11.1%). The calculated risk using the European System for Cardiac Operative Risk Evaluation II was 3.5% (IQR: 2.0% to 8.1%).

Results. Surgery was performed at a median 4 days (IQR: 3 to 6) after NOAC withdrawal. Reduced renal function was predictive for length of intensive care unit stay and administration of red blood cells (p < 0.0001 and p = 0.0291, respectively). The NOAC withdrawal interval significantly influenced postoperative drainage volume

withdrawal of 10 days, compared with 4.2 days without termination. Thirty-day mortality was 3.7%.


Conclusions. A lengthy NOAC withdrawal period, particularly for patients with reduced renal function, is essential for safe open-heart surgery. We conclude that despite official recommendations, patients should whenever possible not be considered for elective cardiac surgery within 10 days of terminating NOAC treatment.

(Ann Thorac Surg 2018;105:702–8) © 2018 by The Society of Thoracic Surgeons

#### **Evolution of PCI in Germany**




#### Mengenentwicklung Koronarchirurgie versus PCI – 2000 bis 2016




Quelle: Darstellung nach J. Cremer, 2016, auf Grundlage von Daten der DGTHG und aus dem Dt. Herzbericht 2000-2014

#### **Evolution of antiplatelet therapy after PCI**





### **2017 ESC Guidelines for STEMI patients**








#### Coronary artery bypass grafting-related bleeding complications in patients treated with ticagrelor or clopidogrel: a nationwide study

Emma C. Hansson<sup>1</sup>, Lena Jidéus<sup>2</sup>, Bengt Åberg<sup>3</sup>, Henrik Bjursten<sup>4</sup>, Mats Dreifaldt<sup>5</sup>, Anders Holmgren<sup>6</sup>, Torbjörn Ivert<sup>7</sup>, Shahab Nozohoor<sup>4</sup>, Mikael Barbu<sup>3</sup>, Rolf Svedjeholm<sup>8</sup>, and Anders Jeppsson<sup>1,9\*</sup>

38% major bleeding!





#### **Ticagrelor Removal From Human Blood**

George O. Angheloiu, MD, <sup>a,b,c</sup> Gabriel B. Gugiu, PнD, <sup>d</sup> Cristian Ruse, PнD, <sup>e</sup> Rishikesh Pandey, PнD, <sup>a</sup> Ramachandra R. Dasari, PнD, <sup>a</sup> Carl Whatling, PнD<sup>f</sup>



#### **VISUAL ABSTRACT**





Angheloiu, G.O. et al. J Am Coll Cardiol Basic Trans Science. 2017;2(2):135-45.

#### HIGHLIGHTS

- Ticagrelor is reversibly bound to albumin.
- CytoSorb and Porapak Q 50-80 mesh remove ticagrelor from bovine serum albumin solution with >99% efficiency.
- CytoSorb removes ticagrelor from human blood and human plasma with >99% efficiency.

...after 10 hours, and 94% after 3-4 hours recirculation

#### **Hypothesis**



- Platelet transfusion is not a definitive solution because of circulating Ticagrelor binding to transfused platelets.
- Using adsorber technology may reduce bleeding complications in patients treated with Ticagrelor.
- Similar effects may be expected in patients treated with NOACs


#### Study design (n = 55)



- Single centre prospective cohort study
- ➤ Patients who underwent emergency cardiac surgery at our institution between June 2016 and June 2018 with preoperative treatment of ticagrelor (n = 43) or rivaroxaban (n = 12).
- Since April 2017 (JACC paper) we routinely installed standardized Cytosorb adsorption into the extracorporeal circulation.


## **Intraoperative setting**





#### Single center prospective cohort study





## **Preoperative data I**



	Cytosorb	Cytosorb	Control	Control	
	Ticagrelor	Rivaroxaban	Ticagrelor	Rivaroxaban	
Demography	(n=32)	<u>(n=7)</u>	(n=11)	(n=5)	p value
age (y)	66	NOSON)	69	72	.33
female (%)	19 🧯	57	18	20	.22
BMI (kg/m²)	27	26	27	27	.79
NYHA class (%)					
	<i>5</i> 6	57	36	60	
III	38	43	64	40	
IV	6	0	0	0	
Comorbidities (%)					
hypertension	91	86	100	100	.99
periph. vasc. disease	28	29	18	40	.68
COLD	38	57	46	60	.92
Renal impairment					.95
normal	28	0	27	0	
moderate	41	71	<i>5</i> 5	80	
severe	31	29	18	20	

## **Preoperative data II**



	Cytosorb	Cytosorb	Control	Control	
	Ticagrelor	Rivaroxaban	Ticagrelor I	Rivaroxaban	
	(n=32)	(n=7)	(n=11)	(n=5)	p value
LVEF (%)		REGAIN A CONTROL			.73
good (>50%)	47	<b>5</b> 7	36	40	
moderate (31-50%)	47	14	64	40	
poor (<30%)	6	29	0	20	
Pathology (%)					
coronary artery dis.	100	100	100	100	
aortic valve disease	9	-	-	-	
mitral valve disease	3	-	9	-	
aortic dissection	3	-	-	-	
atrial fibrillation	13	100	18	100	
EuroSCORE II (%)	3.1	3.9	3.1	3.3	.56
Emergency (%)	100	100	100	100	

#### **Results**

	Cytosorb Ticagrelor	Cytosorb Rivaroxaban	Control Ticagrelor	Control Rivaroxaban	
	(n=32)	(n=7)	(n=11)	(n=5)	p value
Surgical procedure		STOSOTO :			
CABG	84	<b>?</b> 100	91	100	
CABG + AVR	9				
CABG + MVR	3		9		
Aortic replacement	3				
Concomitant surgery					
Afib ablation	9	43	9	40	
LAA occlusion	6	14	9	40	
Time-related outcome					
CPB time	115	80	108	97	.41
X clamp time	77	81	64	70	.54
Total duration	288	184	353	309	.004


### **Bleeding / Length of stay**



	Cytosorb Ticagrelor (n=32)	Cytosorb Rivaroxaban (n=7)	Control Ticagrelor (n=11)	Control <mark>Rivaroxaban</mark> (n=5)	p value
Rethoracotomy	0	0	36	40	.0003
Drainage volume (24hrs)	350	390	890	600	.004
Days in ICU	2	2	3	6	.01
Total length of stay (days)	11	11	14	18	.02

#### Postoperative drainage volume






p = 0.004

#### **Summary**





**Conclusions:** Cytosorb adsorbtion is a safe and effective method to reduce bleeding complications during emergency open-heart surgery in patients with Ticagrelor or Rivaroxaban medication.

#### **Conclusions**



# Both medical and economic benefits of using Cytosorb in Ticagrelor- and Rivaroxaban-loaded patients:

- > reduced operation time
- decreased use of blood products
- saves costs by faster discharge of patients from ICU

The data show that the strategy is a safe and effective method to

- reduce bleeding complications
- improve surgical outcome significantly.

#### The story continues...



	Cytosorb	Cytosorb
	Ticagrelor	<b>Ticagrelor</b>
	(n=32)	(n=61)
Rethoracotomy		
	0	1 (1.6%)
Drainage volume/ 24hrs		
<del>-</del>	420 ± 246	487 ± 222
Days in intensive care		
	3 ± 2	2 ± 3
Total length of stay		
	12 ± 7	12 ± 5
30-days-mortality, n (%)		
	0	1 (1.6%)
Transfusion of platelets		
	11 (34.4%)	20 (32.8%)
Transfusion of red blood cells		
	7 (21.9%)	13 (21.3%)

From 01/2020-05/2020 cytosorb adsorption in 18% of all pats.

#### **Alternatives: Ticagrelor / NOAC antidots**



**PB2452** 

(PhaseBio Pharma. Inc.)

for Ticagrelor

monoclonal ab, phase I (2019)

**Idarucizumab** (Praxbind)

for Dabigatran

5g (100ml)

US\$ 8,385.20

**Andexanet alfa (Ondexxya)** 

for Apixaban, Rivaroxaban

US\$ 2,873.38/ 100 mg vial

low dose: 400 mg bolus + 480 mg iv

total US\$ 25,850.90

US\$ 5,744.38/ 200 mg vial

high dose: 800 mg bolus + 960 mg iv

total US\$ 48,828.42

### **Clinical perspective**



**Standard procedure:** 



Cytosorb adsorption in all ticagrelor and/or

**NOAC-loaded patients during emergency cardiac surgery** 



# **NOACs and CytoSorb**

C. Michael Gibson, M.S., M.D.

Interventional Cardiologist

Professor of Medicine Harvard Medical School

President & CEO of Non-Profit Baim Institute

Founder, Editor-In-Chief www.wikidoc.org





#### **Disclosure**

- Dr. Gibson has received research grant support and consulting fees in the past from all major manufacturers of antiplatelets and antithrombins
- This is an educational lecture and is not intended to be an inducement to use any drug or drug in a fashion that is inconsistent with the drug or device label. Rivaroxaban is not approved for use in acute coronary syndromes in the US, but is so in many other countries
- The slides were prepared by C. Michael Gibson, M.S., M.D. and / or were under the editorial control of C. Michael Gibson, M.S., M.D.

#### **Disclosures**

#### **Present Research/Grant Funding**

Angel Medical Corporation Janssen/Johnson & Johnson CSL Behring SCAD Alliance Baim Institute

Patents and Stocks: None

Equity: nference, Inc.

#### Consultant

Amarin

**Angel Medical Corporation** 

AstraZeneca

Bayer

**Boston Clinical Research Institute** 

**Caladrius Biosciences** 

Cardiovascular Research Foundation

CeleCor Therapeutics

Eli Lilly

**Eidos Therapeutics** 

Gentech

Janssen/ J&J

Kiniksa Pharmaceuticals

MD Magazine

Medtelligence

The Medicines Company

Micodrop, LLC

Microport

**MJHealth** 

Novo Nordisk

Pfizer

**Somahlution** 

**Thrombolytic Science** 

**Verseon Corporation** 

WedMD

#### Consultant

(with monies paid to hospital)
Bayer Corporation
Janssen Pharmaceuticals

#### **Spouse:** Employee of Boston Clinical Research Institute in which she has equity position

Amarin
Amgen
AstraZeneca
Bayer/Janssen/ J&J
Boehringer Ingelheim
Boston Scientific
Cardiovascular Research

Foundation

Calardius Biosciences CeleCor Therapeutics

Chiesi

CSL Behring

DCRI Eidos Therapeutics

Eli Lilly GE Healthcare

Gilead Sciences. Inc.

Gilead Sciences, Inc.

Impact Bio, LTD Kiniksa Pharmaceuticals

MD Magazina

MD Magazine

The Medicines Company

Medtelligence
Medlmmune
The Medicine's Co.
Merck & Co. Inc.
Micodrop, LLC
Microport
MJHealth
Novo Nordisk
PERT Consortium

Pfizer
PharmaMar
Samsung
Sanofi
Somahlution
Syneos Health
Thrombolytic Science

**Verseon Corporation** 

Check the label in your country. Rivaroxaban is not FDA approved in the ACS setting or in patients with atrial fibrillation undergoing stent placement. It is in many other countries. Check your local label. The use of Rivaroxaban in chronic CAD is under regulatory review and is off label at present.


Slide by C. Michael Gibson, M.S., M.D.

#### **Indications for NOAC**

- Atrial fibrillation
- Surgical VTE prophylaxis (knee and hip surgery)
- Medically ill VTE prophylaxis
- VTE treatment (DVT and PE)
- Acute coronary syndrome (Ex US)

#### **Atrial Fibrillation Prevalence Continues to Grow**


Atrial fibrillation prevalence is substantial and expected to grow.




Slides by C. Michael Gibson

#### **Atrial Fibrillation Prevalence Increases with Age**


Atrial fibrillation prevalence increases with age.



# **NOAC** in Atrial Fibrillation



# NOAC in Atrial Fibrillation Major Bleeding



# NOAC Are the Standard of Care For Venous Thromboembolism and Atrial Fibrillation











# Anticoagulants Market to be Worth US\$ 40,158.4 Million by 2026, Says TMR



NEWS PROVIDED BY

Transparency Market Research → May 03, 2019, 06:00 ET SHARE THIS ARTICLE











ALBANY, New York, May 3, 2019 /PRNewswire/ -- TMR's analysts estimate that the global anticoagulants market is expected to touch US\$ 40,158.4 mn by the end of the forecast period. The market was valued US\$ 21,759.3 mn in 2018. The growth of the market is anticipated to occur at a promising 8.0% CAGR during 2018-2026.

From the perspective of drug class, factor Xa inhibitors segment is gaining traction in the global anticoagulants market due to its high usage in various indications such as stroke, heart attack, pulmonary embolism (PE), angina, surgery, and deep venous thrombosis (DVT). On the regional front, North America showcases the highest share in the global anticoagulants market with growing number of several surgical procedures such as knee and hip replacements, and rising healthcare expenditures.

<b>ED Visits for ADEs</b>
By Drug Class (2005-2014)
ED Visits National Estimate, % (95% CI)

2009-2010

18.6 (16.7-20.6)

11.2 (8.2-14.1)

2.5 (1.2-3.9)

4.6 (2.9-6.2)

3.0 (2.6-3.4)

12.0 (9.1-14.9)

3.2 (2.8-3.6)

7.2 (6.6-7.8)

2.9 (2.2-3.7)

3.6 (3.1-4.2)

2011-2012

17.9 (16.6-19.3)

13.3 (11.1-15.4)

2.4 (1.3-3.4)

4.8 (3.2-6.4)

3.1 (2.5-3.7)

12.0 (9.1-14.8)

3.2 (2.7-3.7)

7.9 (7.2-8.5)

3.2 (2.4-4.0)

3.4 (2.8-3.9)

2013-2014

16.1 (14.4-17.8)

17.6 (14.2-21.0)

3.0 (1.6-4.3)

6.6 (4.7-8.5)

2.7 (2.1-3.2)

13.3 (10.8-15.8)

2.8 (2.4-3.2)

6.8 (6.3-7.4)

3.5 (2.6-4.4)

3.0 (2.4-3.5)

JAMA 316.20 (2016): 2115-2125

<b>ED Visits for ADEs</b>
By Drug Class (2005-2014)

2007-2008

18.6 (16.8-20.4)

9.9 (7.6-12.2)

1.6 (0.8-2.4)

4.6 (2.7-6.5)

3.0 (2.5-3.5)

12.8 (9.1-16.6)

3.5 (3.0-3.9)

7.9 (7.1-8.8)

2.5 (1.9-3.0)

3.2 (2.8-3.6)

ED Visits for ADEs	
By Drug Class (2005-2014)	
	Ą

2005-2006

19.2 (17.7-20.7)

7.3 (4.6-9.9)

1.8 (0.8-2.8)

3.7 (2.1-5.2)

2.8 (2.2-3.4)

10.9 (7.3-14.5)

4.1 (3.4-4.8)

7.7 (6.9-8.6)

2.4 (1.9-2.9)

3.2 (2.7-3.7)

**Drug Class** 

**Antibiotics** 

**Anticoagulants** 

Antiplatelets

**Antipsychotics** 

Diabetes agents

Opioid analgesics

Sedative/hypnotic agents


Slide by C. Michael Gibson, M.S., M.D.

**RAS** inhibitors

**NSAIDs** 

Antineoplastic agents

## Increased Morbidity and Mortality in Patients Admitted on a **NOAC** As Compared to Non-Anticoagulated Patients



<sup>1.</sup> Truven, MarketScan Commercial, Medicare Supplemental, last 12 months ending April 30,2015. Medicaid accounts for ~5% of the total bleed related admissions.
2. LOS = The LOS in the Truven report varies by payor. In the YTD 10/2014 report the LOS were 8.0 (12.0), 7.1(9.3) for Commercial, Medicare respectively.
3. The data for mortality from major bleeds ranges from 5.1% (DRESDEN Registry) to 33% (RIETE Registry). Other data such as the ARISTOTLE trial (Granger et al, NEJM 2011) suggest 11–15%.
4. Amin et al. J Manag Care Spec Pharm. 2015;21(10):965-72

# **Management for Patients on NOAC Today**

NOAC	Treatment for Life-Threatening Bleeding
Pradaxa® (dabigatran)	<b>Praxbind</b> (idarucizumab)
<b>Xarelto</b> ® (rivaroxaban)	
Eliquis® (apixaban)	Andexxa (andexanet alfa)
<b>Savaysa</b> ® (edoxaban)	

<sup>1.</sup> Praxbind® prescribing information. Ridgefield, CT. Boehringer Ingelheim Pharmaceuticals, Inc. 2015.
2. Mo Y, Yam FK. Recent advances in the development of specific antidotes for target-specific oral anticoagulants. *Pharmacotherapy*. 2015;35:198-207.
3. Cytosorb. Instructions for use.

<sup>5.</sup> AndexXa™ prescribing information. South San Francisco, CA. Portola Pharmaceuticals, 2016.

### **Andexanet Limitations**

- Approved for reversal of ongoing life threatening bleeding
- NOT approved for reversal in a patient with no bleeding who is to undergo surgery
- Numeric excess of thrombotic events following reversal
- Cost

# **Management for Patients on NOAC in the Future**


NOAC	Prevention of Bleeding	Treatment for Life-Threatening Bleeding
Pradaxa® (dabigatran)		<b>Praxbind</b> (idarucizumab)
<b>Xarelto</b> ® (rivaroxaban)	CytoSorb <sup>®</sup>	
<b>Eliquis</b> ® (apixaban)	Cytosorb	Andexxa (andexanet alfa)
<b>Savaysa</b> ® (edoxaban)		

# **Preventing Bleeding in Cardiac Surgery**








# Cardiac Surgery vs. Non-Cardiac Surgery



# In-vitro Removal of Rivaroxaban by CytoSorb








"Within 1 hour 91.6% of circulating rivaroxaban was removed."

"For normal therapeutic concentrations below 300 µg/L, we expect plasma concentration to be reduced below the critical threshold in 30 to 60 minutes."

# Ticagrelor + Rivaroxaban Removal off-pump

58 y/o male, high bleeding risk, undergoing urgent OPCAB

- PCI dissected LAD. On aspirin, ticagrelor, rivaroxaban (Afib)
- Ongoing chest pain, (+) Tnl
- Urgent OPCAB recommended
- Cytosorb started 1 hour prior to surgery and continued for 1.5 h into CABG
- Cytosorb integrated in hemoperfusion mode
- Operative course uneventful without excess bleeding
- Patient well at 6 months f/u



Dual antithrombotic removal (TIC + RIV) without CPB

# In-vitro Removal of Dabigatran by CytoSorb




Contents lists available at ScienceDirect

International Journal of Cardiology

journal homepage: www.elsevier.com/locate/ijcard

Removal of dabigatran using sorbent hemadsorption Alexandra A. Angheloiu A. George O. Angheloiu b.c.\*



BSA-Dabigatran


Sorbent Column

Infusion Pump

"Dabigatran is robustly removed by a sorbent hemadsorption method already proven successful for ticagrelor. Dabigatran removal restores the aPTT, suggesting reversal of the anticoagulant effect of this drug."


**Volume of column** 


# In-vitro Removal of Edoxaban by CytoSorb



Alexandra A. Angheloiu<sup>1</sup> · Yanglan Tan<sup>2,3</sup> · Cristian Ruse<sup>4</sup> · Scott A. Shaffer<sup>2,3</sup> · George O. Angheloiu<sup>5</sup>

and Albumin Solution





"Sorbent-mediated technology may represent a viable pathway for edoxaban removal from human plasma."

# **Hospital-wide Applications**

- Off-pump cardiac surgery
- Cardiac electrophysiology procedures
- Neurosurgical procedures
  - Reduce risk of life-threatening bleeding
  - Avoid surgery (subdural hematomas)
- Acute stroke (NOAC = contraindication for t-PA)
- Urgent orthopedic procedures
- Urgent GI or oncological procedures
- Trauma

# **Summary: NOAC and CytoSorb**

- NOAC are the standard of care for chronic anticoagulation (Afib, VTE, etc.)
- Aging of the population will only increase use
- Patients on NOAC present unique challenges in the hospital setting due to bleeding risk
- CytoSorb<sup>®</sup> is the only strategy that can prevent bleeding in these patients
- Clinical evidence supports its use for on pump cardiac surgery
- Future studies can establish its use throughout the hospital (any OR, ED, etc.)

# **Closing Remarks**

Efthymios N. Deliargyris, MD, FACC, FESC, FSCAl Chief Medical Officer CytoSorbents Corporation

# Regulatory Status & Clinical Activities

### **European Union**

 CytoSorb is CE Mark label approved for ticagrelor and rivaroxaban removal

### USA

- Ticagrelor removal granted
   FDA Breakthrough Designation
- Ongoing discussions with FDA to set regulatory pathway



Instructions For Use CytoSorb® 300 mL Device

### 1. INTRODUCTION

### 1.1. Intended Use

The CytoSorb Device (CytoSorb) is a non-pyrogenic, sterile, single use polymer based adsorption system designed for the dialysis of physiological fluids in the area of extracorporeal therapies.

#### 1.2. Indications

CytoSorb beat comproven to remove the following:

P2Y <sub>12</sub> Inhibitor	Inflammatory Mediators	Other Substances
Ticagrelor	Cytokines	Bilirubin
Factor Xa Inhibitor		Myoglobin
Rivaroxaban		


CytoSorb is indicated for use in conditions where elevated levels of cytokines and/or bilirubin and/or myoglobin exist.

CytoSorb is indicated for use intraoperatively during cardio-pulmonary bypass surgery for the removal of P2Y, -Inhibitor Ticag-relor and/or Factor Xa-Inhibitor Rivaroxaban.

Study	Design	Region	Start
TISORB  Cagrelor Cyto CORE Hemoadsorption	Prospective, open label trial	United Kingdom	Q3 '20 Continued
CYTATION The OytoSorb® Ticagrelor HemoAdsorp TION Study	Prospective, open label trial	Germany	Q3 '20
STAR Safe and Timely Antithrombotic Removal	Prospective, international registry	Phase 1: EU Phase 2: US + ROW	Q4 '20

# Total Addressable Market For Ticagrelor Removal


### **United States**



- CABG: Coronary Artery Bypass Graft surgery
- AA: Aortic Aneurysm Repair



# **US Market - Sequential Growth**





# **Final Thoughts**

- Antithrombotic drug removal with CytoSorb is a novel solution to a very large unmet hospital need
- Currently no available therapies to <u>prevent</u> bleeding
  - Andexxa® or PB2452 (not yet approved) intended only for use <u>after</u> life-threatening bleeding
- CytoSorb antithrombotic removal in cardiac surgery is safe, effective, easy to implement and is expected to lead to substantial cost savings (dominant value proposition)
- Already approved in E.U. for cardiac surgery (ticagrelor + rivaroxaban) and Breakthrough Designation granted by FDA (ticagrelor)
- Ongoing clinical projects to establish removal of additional NOACs and hospital-wide clinical use
- Market opportunity for all "at-risk" surgeries exceeds \$1.5
   Billion annually in the U.S. alone

### THANK YOU FOR YOUR ATTENTION

